

# **Project Setup Steps**

## **1. Set Up the Raspberry Pi 4B**

### **1.1 Install Raspbian OS**

1. Download the latest Raspbian OS image from the [official Raspberry Pi website](<https://www.raspberrypi.org/downloads/raspberry-pi-os/>).
2. Use software like Balena Etcher to flash the OS image onto a microSD card.
3. Insert the microSD card into the Raspberry Pi and power it up.

### **1.2 Initial Configuration**

1. Connect the Raspberry Pi to a monitor, keyboard, and mouse.
2. Complete the initial setup wizard:
  - Configure Wi-Fi
  - Set the locale
  - Update the software

### **1.3 Install VNC**

1. Install VNC Viewer on your system using browser to start the OS
2. Using Raspberry pi's IP Address connect raspberry pi to VNC Viewer to start the OS

## **2. Install Necessary Software and Libraries**

### **2.1 Update the System**

`sudo apt update`

`sudo apt upgrade`

### **2.2 Install Python and Pip**

`sudo apt install python3 python3-pip`

### **2.3 Install OpenCV**

`sudo apt install libopencv-dev python3-opencv`

### **2.4 Install TensorFlow and Keras**

`pip3 install tensorflow keras`

### **2.5 Install Additional Libraries**

`pip3 install numpy scipy imutils`

`sudo apt install espeak`

### **3. Set Up the Camera**

#### **3.1 Enable the Camera Interface**

1. Open the Raspberry Pi configuration tool:

```
sudo raspi-config  
libcamera-hello
```

2. Navigate to `Interfacing Options` and enable the camera.

3. Reboot the Raspberry Pi:

```
sudo reboot
```

### **4. Connect the Ultrasonic Sensor**

#### **4.1 Wiring**

1. Connect the VCC pin of the ultrasonic sensor to the 5V pin of the Raspberry Pi.
2. Connect the GND pin to a ground pin on the Raspberry Pi.
3. Connect the TRIG pin to a GPIO pin (e.g., GPIO 23).
4. Connect the ECHO pin to another GPIO pin (e.g., GPIO 24).

#### **4.2 Install GPIO Library**

```
sudo apt install python3-rpi.gpio
```

### **5. Download and Set Up the SUPERVISION Code**

#### **5.1 Download Code File**

1. Download the code manually

[https://drive.google.com/drive/folders/1inbCx8TVoCqqKD\\_S0SUUODfJ1cdpVXyd?usp=sharing](https://drive.google.com/drive/folders/1inbCx8TVoCqqKD_S0SUUODfJ1cdpVXyd?usp=sharing)

#### **5.2 Prepare the Model**

1. Download the pre-trained MobileNetV2 model and COCO names file.
2. Ensure they are placed in the correct directory as specified in your code.

## **6. Run the Script and Bluetooth**

Run the script using `python object_detection1.py --labels labels.txt --model mobilenet_v2.tflite`

Using Bluetooth Option connect the earphones to the Bluetooth

## **7. The User setup**

7.1 Wear the Prototype to the user carefully and connect Bluetooth for voice commands

7.2 The Scripts runs for 1 min

7.3 To rerun the script please run the script command again

## **8. Close the Script**

Close the terminal and shut down the system.